



Exhibit E






A service of the National Library of Medicine
and the National Institutes of Health


www.pubmed.gov


[My NCBI](#)
[\[Sign In\]](#) [\[Regis\]](#)


All Databases
PubMed
Nucleotide
Protein
Genome
Structure
OMIM
PMC
Journals
Book


Search PubMed  for

☒ Limits
 ☐ Preview/Index
 ☐ History
 ☐ Clipboard
 ☐ Details

Display AbstractPlus  Show 20

 Sort by

 Send to

All: 1 Review: 0 

Links

☐ 1: [Prog Clin Biol Res. 1987;231:443-57.](#)

Recovery of endotoxin from human plasma by acid oxidative treatments as monitored by an automated microtiter plate-chromogenic substrate Limulus amoebocyte lysate (LAL) assay method.

Tsuji K, Martin PA, Gaunnac GL.

Methods to detect endotoxin in human plasma for LAL testing have been developed by use of a perchloric acid (PCA) or a trifluoroacetic acid (TFA) treatment. These acid oxidative treatments eliminate interferences present in plasma. The mode of action of the acids for the elimination of the interferences is indicated to be by their oxidative actions. The treated plasma is assayed by use of an automated microtiter plate-chromogenic substrate LAL assay method employing a modified Cetus Pro/Pette EXPRESS system. The automated system assays 96 samples in approximately 22 minutes. Assay sensitivity of the automated Pro/Pette system is approximately 0.01 endotoxin units (EU) per ml with a linear range from 0.03 to 0.14 EU/ml. The relative standard deviation (RSD) of the assay system is approximately 6%. The improved precision of the LAL assay by the automated system has been attributed to its fast speed (2 minutes) in dispensing the unstable LAL reagent to all 96 samples in one assay block. The total operation, including the acid oxidative treatments and the 96 sample LAL assay, takes approximately 45 minutes. The recovery of spiked endotoxin from plasma ranged from 75 to approximately 100% by adding 0.4 ml of either 0.1-0.2 N PCA or 0.2-0.35 N TFA to 0.2 ml plasma and reacting for 15 minutes at 37 degrees C. The standard deviation of the assay is approximately 0.1 EU/ml plasma. The method has been proven to be rugged, simple, rapid, and cost effective thus would be suitable for clinical application.

PMID: 3035582 [PubMed - indexed for MEDLINE]

Related Links

Automation of chromogenic substrate Limulus amoebocyte lysate assay method for endotoxin by robotic system [Appl Environ Microbiol. 1984]

Evaluation of robot automated chromogenic substrate LAL endotoxin assay method for pharmaceutical production testing [J Pharm Biotechnol. 1985]

Automation and computerization of chromogenic LAL assay method for bacterial endotoxin using 96-well microtiter plates [J Automat Sci Technol. 1990]

Endotoxin testing in blood. [Prog Clin Biol Res. 1987]

Optimization of detection of bacterial endotoxin in plasma with the Limulus test. [J Lab Clin Med. 1990]

See all Related Articles...

Display AbstractPlus  Show 20

 Sort by

 Send to

[Write to the Help Desk](#)

[NCBI](#) | [NLM](#) | [NIH](#)

[Department of Health & Human Services](#)

[Privacy Statement](#) | [Freedom of Information Act](#) | [Disclaimer](#)